

Lateral "J"

Ideal Application

The tubular lateral "J" buttress is ideal for supporting lateral patellar subluxation, dislocations or lateral patellofemoral malalignment.

Features and Benefits

- Tubular lateral "J" buttress for patellar stabilization
- Lateral to medial strapping for additional stabilization
- 1/8" neoprene for heat and compression or Drytex material for patients with allergies or irritations to neoprene
- Available in left or right configuration

Ordering Information

PART NUMBER	
11-0320-x-xxxxx	Right
11-0321-x-xxxxx	Left
11-0659-x-06060	Drytex® Right
11-0660-x-06060	Drytex® Left
<i>SIZES AVAILABLE: XS-XXL Measurements taken 6" above mid-patella</i>	
x=1=XS=13"-15 1/2" (33-39cm)	x=2=S=15 1/2"-18 1/2" (39-47cm)
x=3=M=18 1/2"-21" (47-53cm)	x=4=L=21"-23 1/2" (53-60cm)
x=5=XL=23 1/2"-26 1/2" (60-67cm)	x=6=XXL=26 1/2"-29 1/2" (67-75cm)
Colors Available: black	

Suggested HCFA Code: L1800 + L2795



NEVER STOP GETTING BETTER™

www.djortho.com